

见证更多欢笑  
Smile  
See More Happiness

CE 1639

# Oxygen Concentrator

## Model:KSOC-10



- ◆ America PSA technology offer nature oxygen
- ◆ Reliable and durable oil free compressor
- ◆ Available 24 hours continuous working
- ◆ Self diagnostic system with error code indication
- ◆ Optonal dual flow and flow splitter for multi people use

### Functions:

Power off alarm, Overload protection, High/Low pressure alarm, Temperature alarm, Error code indication, Optional nebulizer,oxygen purity alarm and Purity display.

### Specification:

Model	Oxygen Purity	Flow Range	Noise	Input Voltage	Output Pressure	Produce Size(MM)	Power	N.W.	G.W.	Paking Size(MM)
KSOC-10	93±3%	0-10L	≤60dB	220V/110V	30-80kPa	380*320*620	820VA	23kgs	26kgs	465*395*700